

INTRODUCTION

- Reversed-phase chromatography (RPC) is popular for peptide purification
- Establishing a chromatographic process is labor- and time-intensive
- Calibrating a mechanistic model can be done with fewer experiments
- Such a model can then be used to quantify the impact of Process Parameters on Quality Attributes and Process Performances
- This provides the basis to determine the Proven Acceptable Range of process parameters

MATERIALS AND METHODS

Experimental details

This case study extends previous work done for RPC using a Kromasil 100A 10 µm C8 medium packed in a 25 x 0.46 cm column (see Guetaz et al. 2013 Journal of Chromatography A).

Mechanistic model

The model uses the mixing cell approach with the Linear Driving Force approximation together with a specialized isotherm for reversed-phase binding.

$$\bar{C}_i = \underbrace{\frac{H_{1,i} C_i}{1 + \frac{H_{1,i}}{\bar{N}_{1,i}} C_i}}_{\text{Langmuir type interaction}} + \underbrace{\frac{H_{2,i} C_i \left(1 + I_{ii} \frac{H_{2,i}}{2\bar{N}_{2,i}} C_i\right)}{1 + \frac{H_{2,i}}{\bar{N}_{2,i}} C_i + I_{ii} \frac{H_{2,i}}{(2\bar{N}_{2,i})^2} C_i^2}}_{\text{Moreau type interaction}}$$

where $\bar{C}_i$ is the concentration in the solid phase, C_i is the concentration in the liquid phase, $H_{1,i}$ and $H_{2,i}$ are the Henry coefficients for Langmuir and Moreau interactions, $\bar{N}_{1,i}$ and $\bar{N}_{2,i}$ are the capacities for the Langmuir and Moreau interactions, and I_{ii} is an adsorbate-adsorbate interaction parameter

RESULTS

Excellent agreement between experiments and simulation for the product

- Experimental and simulated data show:
 - Langmuir type, tailing elution peak for low loading amounts (≤ 0.2 g/L column)
 - Moreau type, fronting elution peak for higher loading amounts




Fig. 1: Prediction of isocratic (25.58 (v/v)% acetonitrile in water) peptide elution over a range of loading amounts. Scatter points show experimental data, lines show simulation predictions.

Excellent agreement between experiments and simulation for competitive binding with impurities




Fig. 2: Prediction of isocratic (25.58 (v/v)% acetonitrile in water) peptide elution for a 7.5 g/L load including impurities. Scatter points show experimental data, lines show simulation predictions.

RESULTS

The following process is considered for preparative purposes:

- 2 Bed Volumes equilibration with 22 v/v% acetonitrile in water
- Variable loading over a range of 1.6 to 6.5 g/L column
- Elution gradient from 22 v/v% up to 36 v/v% over a range of 40-160 minutes

Quantification of the impact of process parameters

- Simulations were performed to explore the Design Space by varying load (in g/L column) and gradient duration (in minutes)
- The collection strategy was defined
 - To: Maximize yield
 - While: Maintaining purity > 97%




Fig. 3: Simulated impact of process parameters.

Recommendation of the Proven Acceptable Range

- Simulations were performed to determine
 - the Edge of Failure (quality constraints: L1<1.1%, L2<0.032%, H1<1.4 %)
 - economically viable conditions (performance constraints: Yield>86%, Productivity>0.03 g/L/min)
- The following Proven Acceptable Range was recommended
 - Load: 2.3 – 5.8 g/L column
 - Gradient duration: 51 – 114 min




Fig. 4: Determination of the Edge of Failure (left) and economically viable conditions (right). Recommendation of the standard operating point (SOP) and Proven Acceptable Range (PAR).

CONCLUSION

- A mechanistic model for reversed-phase peptide chromatography was developed. The model was proven capable of representing complex peptide interactions with the solid phase, leading to peculiar behaviors.
- The model was then used to identify the range of process parameters allowing to meet both quality and performance constraints.
- This provides a solid basis to determine the Edge of Failure and provide recommendations on the Proven Acceptable Range.

VISIT OUR WEBSITE

